

Evaluated By: \_\_\_\_\_  
 Mark For: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Account # : \_\_\_\_\_

Company: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Phone: \_\_\_\_\_

**AVENGER SCOOTERS - Please Choose One**

E1230

GA-531 **Avenger - 3 wheel**  
 3 Wheel - Weight Capacity: 500 lbs.  
**\$3,495.00**

GA-541 **Avenger - 4 wheel**  
 4 Wheel - Weight Capacity: 500 lbs.  
**\$3,595.00**

**COLOR - Please Choose One**

Silver N/C

**SEAT COLOR - Please Choose One**

Gray Vinyl N/C

**STANDARD SEATING - Match with Corresponding Scooter**

22" Wide X 17.5" Deep N/C

**STANDARD TIRES - Please Choose One**

GA-531 - **Avenger** - 3 wheel  
 (1) 10" Pneumatic Front Tire N/C  
 (2) 13" Pneumatic Rear Tires N/C

GA-541 - **Avenger** - 4 wheel  
 (2) 13" Pneumatic Front Tires N/C  
 (2) 13" Pneumatic Rear Tires N/C

**COMPLETE SET OF FOAM FILLED FLAT FREE TIRES - Additional Charges Apply**

GA-531 Avenger - 3 wheel **\$375.00**

GA-541 Avenger - 4 Wheel **\$500.00**

**BATTERIES (2 Needed to Operate Companion)**

(2) 22NF (55 Amp) **\$260.00 pair**

**STANDARD FEATURES**

Battery Charger: On Board - 4 Amp  
 Flip Up & Adjustable Arm Rests  
 Infinite Adjustable Tiller  
 Front Basket  
 Safety Headlights & Tail lights  
 Swivel Seat (360 Degrees)  
 Self Diagnostics  
 Sleep Mode

**OPTIONAL ACCESSORIES - Additional Charges Apply**

Flag **\$16.00**     Stump Support **\$226.00**     Quad Cane Holder **\$97.00**     Crutch Holder Tube **\$70.00**  
 Seat Belt **\$38.00**     Pack n' Go **\$52.00**     Rear View Mirror **\$38.00**     Accessory Bracket\* **\$32.00**  
 Cup Holder **\$35.00**     Scooter Cover **\$98.00**     Cane Holder Tube **\$40.00**     Oxygen Tank Holder\* **\$110.00**

**SPECIFICATIONS**

Model	Length	Width	Turning Radius	Speed	Range (miles)
Avenger - 3 wheel	51.5"	24.63"	42"	7 mph	up to 20
Avenger - 4 wheel	57.5"	24.63"	48.5"	7 mph	up to 18

**WARRANTY**

**Limited Lifetime**    **Two Year**  
 Metal Frame    Drive Train  
 Battery Charger

ASSIGNMENT OF A HCPCS CODE DOES NOT INDICATE COVERAGE OR PAYMENT FOR THE ITEM.

Golden Technologies reserves the right to alter any specification or MSRP without prior notice!

Avenger092806

**Totals**

\$

\$

\$

\$

\$

\$

**Total Cost**