

Softouch  
with STAIN DEFENSE

Evergreen



Cabernet



Admiral



Palomino



Dorado Tan



C O M F O R T E R S E R I E S

# Rise



## Tall

Featuring our extra soft seat and back pillows with lumbar support and our "grip rite" arm design.

Recommended for persons 6'3" to 6'8".



Rise.



Roll.



Rest.



Ride.



**GOLDEN**  
technologies

YOUR LIFE IN MOTION

# Tall

## C O M F O R T E R S E R I E S

Model Number	PR-501T
Size	TALL
Weight Capacity	375 LB
Recline Positions	3
Standard Fabrics	5
Fabric Options	32
Chaise Seat	YES
Back Type	SEAM
Head Rest Cover	YES
Arm Covers	YES
Battery Back Up	YES
Overall Width	34.5"
Overall Height	47"
Width Between Arms	23"
Floor to Top of Seat	22.5"
Seat Depth	24"
Seat to Top of Back	30"
Seat to Top of Arm	6"
Floor to Back Edge of Seat Raised	31"
Floor to Front Edge of Seat Raised	26"
Seat Angle Raised	18°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	36"
Back Angle Closed	115°
Back Angle Reclined	145°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	19"
Extension Reclined from Back of Seat to End of Foot Rest	45"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YRS.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YRS.
Warranty on Electrical Labor	1 YR.
Removable Back	YES
Customizable	YES
Left Hand Control Available	YES
Extra Pocket Available	YES
Smart Tek System	YES

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

### Safety System

Batteries not included. Safety System is designed to bring chair from recline to seated position in the event of a power outage.

Hand-blown, Overstuffed, Conjugated Microfiber Dacron Fill  
**SUPER SUPPORTIVE**

Northern Maple Wood  
• Double-doweled at Stress Points  
• Reinforced Miter Blocks  
• Stapled and Glued  
**STRONGEST IN THE INDUSTRY**

2" of Block Foam  
**EXTRA COMFORT**

Quick Release Hand Control  
**EASY DISCONNECT**


Rear Landing Legs  
**EXTRA STABILITY**

5 Heavy Gauge Springs  
**LONG LASTING FOR INCREDIBLE COMFORT**

Widest Fabric Selection in Industry  
**FITS IN ANY HOME DÉCOR**

10 Gauge Steel Scissor Mechanism  
**STRONG & LONG LASTING**

Adjustable Anti-Skid/Scrape Feet  
**STABILITY & SUPPORT**



501 M shown



501 M shown



501 M shown



501 M shown

### Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.

### Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.

### Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



### Smart Tek

Our luxury lift & recliners are equipped with electronics able to assure a secure lift every time you need it.



### Custom Sizing...

One size doesn't fit all with Golden Technologies seat lift chairs.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Our chairs meet some of the toughest safety standards ever. Our entire chair meets ETL (Electrical Testing Lab) standards to ensure electrical safety. ETL is a nationally recognized, third party testing facility.

