



CAMBRIDGE MEDIUM/LARGE

\$1,349.00

The Cambridge power lift recliner features a supportive biscuit backrest design as well as a soft, comfortable seat for hours of relaxation. [It's also available in Small/Medium!](#)

SKU: PR401-MLA

Categories: [1 Zone Chairs](#), [Heavy Duty](#), [Lift Recliners](#), [Signature](#)

ADDITIONAL INFORMATION

Size	Med/Large
Class II Medial Device	Yes
Weight Capacity	375 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	Standard

Overall Width	35"
Overall Height	44.5"
Seat Width	22"
Floor to Top of Seat	22"
Seat Depth	22"
Seat to Top of Back	31"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	19"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard

USB Charging Port	Hand Control
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA