



CAMBRIDGE SMALL MEDIUM

\$1,249.00

The Cambridge power lift recliner features a supportive biscuit backrest design as well as a soft, comfortable seat for hours of relaxation. [It's also available Medium/Large!](#)

SKU: PR401-SME

Categories: [1 Zone Chairs](#), [Heavy Duty](#), [Lift Recliners](#), [Signature](#)

ADDITIONAL INFORMATION

Size	Small/Med
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	Standard



Overall Width	33"
Overall Height	41.5"
Seat Width	20"
Floor to Top of Seat	19"
Seat Depth	19"
Seat to Top of Back	27"
Seat to Top of Arm	6.5"
Distance Required from Wall Reclined	18 "
Extension Reclined from Back of Seat to End of Foot Rest	41.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	No
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	No
Straight Lift Option	Standard

USB Charging Port

Hand Control

Battery Backup

9 volt

**SmartTek Diagnostic
System**

Standard

Place of Origin

Imported