

FDA Class II Medical Device



## CAMBRIDGE SMALL

**\$1,599.00**

Enjoy the simplicity of comfort in the Cambridge. Available in Small, Medium & Large. 1-Zone Comfort Zone Positioning System.

---

**SKU:** PR401-SMA

**Categories:** [1 Zone Chairs](#), [DeLuna](#), [Lift Recliners](#), [Signature](#)

## ADDITIONAL INFORMATION

<b>Size</b>	Small
<b>Class II Medical Device</b>	Yes
<b>User Weight Capacity</b>	<a href="#">400 lbs</a>
<b>Positioning</b>	3 Position
<b>Comfort Zone(s)</b>	1
<b>Twilight Tilt Technology</b>	No
<b>ZG+</b>	No
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	No
<b>Back Type</b>	Biscuit
<b>Removable Back Type</b>	Quick Release
<b>Head Rest Cover</b>	Standard
<b>Arm Rest Covers</b>	Standard



<b>Overall Width</b>	31"
<b>Overall Height</b>	40.5"
<b>Seat Width</b>	20"
<b>Floor to Top of Seat</b>	18"
<b>Seat Depth</b>	18.5"
<b>Seat to Top of Back</b>	27"
<b>Seat to Top of Arm</b>	7.5"
<b>Distance Required from Wall Reclined</b>	15"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	41"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	No
<b>Heat &amp; Vibrating Massage Option</b>	No
<b>HeatWave Technology Option</b>	No
<b>Articulating Headrest Option</b>	NO
<b>Power Adjustable Lumbar Option</b>	No
<b>Foot Rest Extension Option</b>	No
<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	NO

<b>Battery Backup</b>	9 volt
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	Imported