

FDA Class II Medical Device



CAMBRIDGE

\$1,599.00 - \$1,699.00 Price range:
\$1,599.00 through \$1,699.00

Enjoy the simplicity of comfort in the Cambridge. Available in Small, Medium & Large. 1-Zone Comfort Zone Positioning System.

SKU: PR401

Categories: [1 Zone Chairs](#), [DeLuna](#), [Lift Recliners](#), [Signature](#)

VARIATIONS

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
<p>FDA Class II Medical Device</p> 	PR401-LAR	\$1,699.00			Enjoy the simplicity of comfort in the Cambridge. Available in Small, Medium & Large. 1-Zone Comfort Zone Positioning System.	Large
<p>FDA Class II Medical Device</p> 	PR401-MED	\$1,649.00			Enjoy the simplicity of comfort in the Cambridge. Available in Small, Medium & Large. 1-Zone Comfort Zone Positioning System.	Medium

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
<small>FDA Class II Medical Device</small> 	PR401-SMA	\$1,599.00			Enjoy the simplicity of comfort in the Cambridge. Available in Small, Medium & Large. 1-Zone Comfort Zone Positioning System.	Small

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Large, Medium, Small
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Quick Release
Head Rest Cover	Standard
Arm Rest Covers	Standard
Overall Width	31"
Overall Height	40.5"
Seat Width	20"
Floor to Top of Seat	18"



Seat Depth	18.5"
Seat to Top of Back	27"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	15"
Extension Reclined from Back of Seat to End of Foot Rest	41"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	No
Heat & Vibrating Massage Option	No
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	No
Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	Imported

