

FDA Class II Medical Device



**CLOUD ELITE WITH TWILIGHT & NIRVANA**

**\$3,299.00**

The PR527 Cloud™ Elite blends a classic design with Golden's overstuffed bucket-back. It offers our patented Twilight Tilt Technology & MaxiComfort® Positioning, plus optional Nirvana Swedish Air Massage with FIR Heat. Available in two sizes - Small/Medium & Medium/Large.

**SKU:** PR527

**Categories:** [5 Zone Chairs](#), [Lift Recliners](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

**Tags:** [New](#)

**VARIATIONS**

**Image**

**SKU**

**Price**

**Stock Status** **Stock Quantity**

**Description**

**Size**

FDA Class II Medical Device



PR527-MLA \$3,299.00

The PR527 Cloud™ blends a classic design with Golden's overstuffed bucket-back. It offers our patented Twilight Tilt Technology & MaxiComfort® Positioning, plus optional Nirvana Swedish Air Massage with FIR Heat. Med/Large

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
<small>FDA Class II Medical Device</small> 	PR527-SME	\$3,299.00			The PR527 Cloud™ blends a classic design with Golden's overstuffed bucket-back. It offers our patented Twilight Tilt Technology & MaxiComfort® Positioning, plus optional Nirvana Swedish Air Massage with FIR Heat. Available in two sizes - Small/Medium & Medium/Large.	Small/Med

### ADDITIONAL INFORMATION

<b>Dimensions</b>	N/A
<b>Size</b>	Med/Large, Small/Med
<b>Class II Medical Device</b>	Yes
<b>User Weight Capacity</b>	<a href="#">400 lbs</a>
<b>Positioning</b>	Twilight Tilt
<b>Comfort Zone(s)</b>	5
<b>Twilight Tilt Technology</b>	Yes
<b>ZG+</b>	Yes
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	No
<b>Back Type</b>	Seam

<b>Removable Back Type</b>	Quick Release
<b>Head Rest Cover</b>	No
<b>Arm Rest Covers</b>	No
<b>Overall Width</b>	38.5"
<b>Overall Height</b>	44"
<b>Seat Width</b>	23.5"
<b>Floor to Top of Seat</b>	19.5"
<b>Seat Depth</b>	20.5"
<b>Seat to Top of Back</b>	29"
<b>Seat to Top of Arm</b>	8.5"
<b>Distance Required from Wall Reclined</b>	21"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	42"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Outside Arm
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	Yes
<b>Nirvana Air Massage &amp; Heat Available</b>	Yes
<b>Articulating Headrest Option</b>	Standard

<b>Power Adjustable Lumbar Option</b>	Standard
<b>Foot Rest Extension Option</b>	No
<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	Right Arm
<b>Battery Backup</b>	Lithium-ion
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA