

FDA Class II Medical Device



CLOUD+ LARGE WITH TWILIGHT

\$3,279.00

Relax & recharge your comfort in the Cloud+, standard with Twilight Tilt & available with Nirvana™ Swedish Air Massage with HeatWave™ Full-Body FIR Heat.

SKU: PR525-LAR

Categories: [5 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Size	Large
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	Twilight Tilt
Comfort Zone(s)	5
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No



Arm Rest Covers	No
Overall Width	39"
Overall Height	48"
Seat Width	25"
Floor to Top of Seat	22"
Seat Depth	22"
Seat to Top of Back	32"
Seat to Top of Arm	8"
Distance Required from Wall Reclined	25"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Compartment
Extra Pocket Option	No
Heat & Vibrating Massage Option	Yes
Nirvana with HeatWave Option	Yes
HeatWave Technology Option	Yes
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	Yes
Straight Lift Option	Standard



USB Charging Port	YES
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA