

FDA Class II Medical Device



CLOUD MEDIUM/LARGE


\$5.00

Relax and recharge in comfort and style with the PR517 Cloud Power Lift Recliner - 4-Zone Comfort Zone positioning system!

SKU: PR517-MLA

Categories: [4 Zone Chairs](#), [Lift Recliners](#), [Patented MaxiComfort Positioning](#)

VARIATIONS

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
 <p><small>FDA Class II Medical Device</small></p>	PR517-MLA	\$5.00				Med/Large

ADDITIONAL INFORMATION

Size	Med/Large
Class II Medical Device	Yes
Weight Capacity	400 lbs

Positioning	MaxiComfort®
Comfort Zone(s)	4
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	No
Overall Width	38"
Overall Height	45"
Seat Width	20"
Floor to Top of Seat	20"
Seat Depth	21.5"
Seat to Top of Back	32"
Seat to Top of Arm	8.5"
Distance Required from Wall Reclined	26.5"
Extension Reclined from Back of Seat to End of Foot Rest	46"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket

Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	Yes
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA