



CLOUD MEDIUM LARGE

\$2,399.00

The Cloud™ is an *extremely* comfortable and popular power lift recliner offering our MaxiComfort® Positioning & custom bucket seat design for ultimate comfort and support.

SKU: PR510-MLA

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Size	Med/Large
Class II Medial Device	Yes
Weight Capacity	375 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Bolt
Head Rest Cover	Standard

Arm Rest Covers	No
Overall Width	38"
Overall Height	45"
Seat Width	20"
Floor to Top of Seat	20"
Seat Depth	21.5"
Seat to Top of Back	32"
Seat to Top of Arm	8.5"
Distance Required from Wall Reclined	26.5"
Extension Reclined from Back of Seat to End of Foot Rest	46"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	Yes
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes

Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA