



## CLOUD MEDIUM WIDE

**\$2,799.00**

Your favorite chair just got an upgrade! The Cloud™ is a popular Golden power lift recliner offering our MaxiComfort® Positioning & 600 lbs. weight capacity.

---

**SKU:** PR510-MXW

**Categories:** [2 Zone Chairs](#), [Favorites](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

**Tags:** [Best Seller](#)

## ADDITIONAL INFORMATION

<b>Size</b>	Medium/Wide
<b>Class II Medial Device</b>	Yes
<b>Weight Capacity</b>	<a href="#">600 lbs</a>
<b>Positioning</b>	MaxiComfort®
<b>Comfort Zone(s)</b>	2
<b>Twilight Tilt Technology</b>	No
<b>ZG+</b>	No
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	No
<b>Back Type</b>	Biscuit
<b>Removable Back Type</b>	Quick Release
<b>Head Rest Cover</b>	Standard
<b>Arm Rest Covers</b>	No



<b>Overall Width</b>	45"
<b>Overall Height</b>	48"
<b>Seat Width</b>	29"
<b>Floor to Top of Seat</b>	21"
<b>Seat Depth</b>	21.5"
<b>Seat to Top of Back</b>	32"
<b>Seat to Top of Arm</b>	7.5"
<b>Distance Required from Wall Reclined</b>	21"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	44"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty on Electrical Parts</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	No
<b>Articulating Headrest Option</b>	NO
<b>Power Adjustable Lumbar Option</b>	No
<b>Foot Rest Extension Option</b>	Yes
<b>Straight Lift Option</b>	Standard



**USB Charging Port**

Right Arm

**Battery Backup**

Lithium-ion

**SmartTek Diagnostic  
System**

Standard

**Place of Origin**

USA