



CLOUD MEDIUM WIDE

\$2,899.00

Your favorite chair just got an upgrade! The Cloud™ is a popular Golden power lift recliner offering our MaxiComfort® Positioning & 600 lbs. weight capacity.

SKU: PR510-MXW

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Size	Medium/Wide
Class II Medial Device	Yes
Weight Capacity	600 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Quick Release
Head Rest Cover	Standard
Arm Rest Covers	No

Overall Width	45"
Overall Height	48"
Seat Width	29"
Floor to Top of Seat	21"
Seat Depth	21.5"
Seat to Top of Back	32"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	21"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard

USB Charging Port

Right Arm

Battery Backup

Lithium-ion

SmartTek Diagnostic System

Standard

Place of Origin

USA