



CLOUD+ LARGE

\$3,099.00

Elevate your senses in the Cloud+, available with our new Nirvana™ Swedish Air Massage with HeatWave™ Full-Body Far Infrared Heat Therapy. Available in 3 sizes.

SKU: PR511-LAR

Categories: [4 Zone Chairs](#), [Favorites](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Size	Large
Class II Medial Device	Yes
Weight Capacity	375 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	4
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No

Overall Width	39"
Overall Height	48"
Seat Width	25"
Floor to Top of Seat	22"
Seat Depth	22"
Seat to Top of Back	32"
Seat to Top of Arm	8"
Distance Required from Wall Reclined	25"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Compartment
Extra Pocket Option	No
Heat & Vibrating Massage Option	No
Nirvana with HeatWave Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard

Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	YES
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA