

FDA Class II Medical Device



CLOUD+ MEDIUM

\$2,999.00

Elevate your senses in the Cloud+, available with our new Nirvana™ Swedish Air Massage with HeatWave™ Full-Body Far Infrared Heat Therapy. Available in 4 sizes.

SKU: PR511-MED

Categories: [4 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medical Device	Yes
Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	4
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No

Overall Width	37.5"
Overall Height	46"
Seat Width	23"
Floor to Top of Seat	21"
Seat Depth	21"
Seat to Top of Back	30"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	22"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Compartment
Extra Pocket Option	No
Heat & Vibrating Massage Option	No
Nirvana With HeatWave Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard

Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	YES
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA