



CLOUD SMALL MEDIUM

\$2,589.00

The Cloud™ by Golden is an extremely comfortable and popular power lift recliner offered in two sizes. These 2-Zone chairs feature our patented MaxiComfort Positioning. Add on our optional, NEW full body Heat Wave Technology with Far-Infrared heat therapy for added well being and relaxation!

SKU: PR510-SME

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Zero Gravity](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

| | |
|----------------------------|---------------------|
| Size | Small/Med |
| Weight capacity | 375 |
| Positioning | MaxiComfort |
| Twilight Tilt | No |
| Chaise seat | Yes |
| Pocketed coil seat | No |
| Back type | Biscuit |
| Removable Back Type | Bolt |
| Head rest cover | YES |
| Arm rest covers | No |

| | |
|-----------------------------------|---------------------|
| Overall width | 34" |
| Overall height | 41" |
| Width between arms | 18" |
| Floor to top of seat | 19" |
| Seat depth | 19" |
| Seat to top of back | 27" |
| Seat to top of arm | 7.5" |
| Dis req from wall reclined | 23" |
| Ext reclined from back | 43" |
| Life time warranty frame | Yes |
| Warranty electrical parts | 3 Yrs. |
| Pro rated 4 thru 7 | Yes |
| Warranty mechanical labor | 3 Yrs. |
| Warranty electrical labor | 1 Yr. |
| Hand control routing | Inside Pocket |
| Extra pocket option | Standard |
| Heat massage option | Yes |
| Heatwave technology option | NEW! Yes |
| Adjustable headrest option | No |
| Power adjustable lumbar | No |
| Foot rest ext option | No |
| Straight lift option | Yes |
| Usb charging port | NO |
| Features | HeatWave Technology |
| Place of origin | USA |