



## CLOUD SMALL MEDIUM

**\$2,589.00**

The Cloud by Golden is an extremely comfortable and popular power lift recliner featuring hybrid sizing. Offered in three sizes, these chairs offer our unique, ergonomic Cloud seating system and custom bucket seat design for long-term comfort and support.

---

**SKU:** PR510-SME

**Categories:** [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Zero Gravity](#)

## ADDITIONAL INFORMATION

<b>Size</b>	Small/Med
<b>Class II Medial Device</b>	YES
<b>Weight Capacity</b>	<a href="#">375</a>
<b>Positioning</b>	MaxiComfort
<b>Comfort Zone(s)</b>	2
<b>Twilight Tilt Technology</b>	NO
<b>ZG+</b>	NO
<b>Chaise Seat</b>	YES
<b>Pocketed Coil Seat</b>	NO
<b>Back Type</b>	Biscuit
<b>Removable Back Type</b>	Bolt



<b>Head Rest Cover</b>	Standard
<b>Arm Rest Covers</b>	No
<b>Overall Width</b>	33.5"
<b>Overall Height</b>	41"
<b>Seat Width</b>	18.5"
<b>Floor to Top of Seat</b>	19"
<b>Seat Depth</b>	19.5"
<b>Seat to Top of Back</b>	27"
<b>Seat to Top of Arm</b>	7.5"
<b>Distance Required from Wall Reclined</b>	23"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	43"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	YES
<b>Warranty on Electrical Parts</b>	3 YR
<b>Pro Rated Years 4 Thru 7</b>	YES
<b>Warranty on Mechanical Labor</b>	3 YR
<b>Warranty on Electrical Labor</b>	1 YR
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	YES
<b>HeatWave Technology Option</b>	YES
<b>Articulating Headrest Option</b>	NO
<b>Power Adjustable Lumbar Option</b>	NO



<b>Foot Rest Extension Option</b>	NO
<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	NO
<b>Battery Backup</b>	9 volt
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA