

FDA Class II Medical Device



CLOUD SMALL MEDIUM

\$2,399.00

The Cloud™ is an *extremely* comfortable and popular power lift recliner offering our MaxiComfort® Positioning & custom bucket seat design for ultimate comfort and support.

SKU: PR510-SME

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

ADDITIONAL INFORMATION

Size	Small/Med
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Bolt
Head Rest Cover	Standard



Arm Rest Covers	No
Overall Width	33.5"
Overall Height	41"
Seat Width	18.5"
Floor to Top of Seat	19"
Seat Depth	19.5"
Seat to Top of Back	27"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	23"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	Yes
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	No



Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA