

FDA Class II Medical Device



CLOUD WITH TWILIGHT & OPTIONAL NIRVANA

\$3,299.00

The PR527 Cloud™ blends a classic design with Golden's overstuffed bucket-back. It offers our patented Twilight Tilt Technology & MaxiComfort® Positioning, plus optional Nirvana Swedish Air Massage with FIR Heat. Available in two sizes - Small/Medium & Medium/Large.

SKU: PR527

Categories: [5 Zone Chairs](#), [Lift Recliners](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

Tags: [New](#)

VARIATIONS

Image

SKU

Price

Stock Status Stock Quantity

Description

Size

The PR527 Cloud™ blends a classic design with Golden's overstuffed bucket-back. It offers our patented Twilight Tilt Technology & MaxiComfort® Positioning, plus optional Nirvana Swedish Air Massage with FIR Heat.

Med/Large

PR527-MLA \$3,299.00



FDA Class II Medical Device

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
 FDA Class II Medical Device	PR527-SME	\$3,299.00			The PR527 Cloud™ blends a classic design with Golden's overstuffed bucket-back. It offers our patented Twilight Tilt Technology & MaxiComfort® Positioning, plus optional Nirvana Swedish Air Massage with FIR Heat. Available in two sizes - Small/Medium & Medium/Large.	Small/Med

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Med/Large, Small/Med
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	Twilight Tilt
Comfort Zone(s)	5
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam



Removable Back Type	Quick Release
Head Rest Cover	No
Arm Rest Covers	No
Overall Width	38.5"
Overall Height	44"
Seat Width	23.5"
Floor to Top of Seat	19.5"
Seat Depth	20.5"
Seat to Top of Back	29"
Seat to Top of Arm	8.5"
Distance Required from Wall Reclined	21"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Outside Arm
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	Yes
Nirvana Air Massage & Heat Available	Yes
Articulating Headrest Option	Standard

Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	No
Straight Lift Option	Standard
USB Charging Port	Right Arm
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA