



CLOUD WITH TWILIGHT SMALL/MEDIUM

\$2,920.00

Experience true weightlessness. The Cloud™ is a luxury power lift recliner offering our patented Twilight Tilt Technology & MaxiComfort® Positioning for long-term comfort and support.

SKU: PR515-SME

Categories: [5 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort Positioning with Tilt](#), [Twilight](#), [Zero Gravity \(MaxiComfort\)](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Size	Small/Med
Class II Medial Device	YES
Weight Capacity	375
Positioning	Twilight
Comfort Zone(s)	5
Twilight Tilt Technology	YES
ZG+	YES
Chaise Seat	YES
Pocketed Coil Seat	NO
Back Type	Biscuit
Removable Back Type	Quick Release



Head Rest Cover	Standard
Arm Rest Covers	No
Overall Width	35.5"
Overall Height	44"
Seat Width	23.5"
Floor to Top of Seat	20"
Seat Depth	18.5"
Seat to Top of Back	28"
Seat to Top of Arm	8"
Distance Required from Wall Reclined	23"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	YES
HeatWave Technology Option	YES
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard

Foot Rest Extension Option	YES
Straight Lift Option	Standard
USB Charging Port	Control Box
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA