

FDA Class II Medical Device



CLOUD WITH TWILIGHT

\$3,439.00

Experience true weightlessness. The Cloud™ is a luxury power lift recliner offering our patented Twilight Tilt Technology & MaxiComfort® Positioning for long-term comfort and support.

SKU: PR515

Categories: [5 Zone Chairs](#), [Lift Recliners](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

Tags: [Best Seller](#)

VARIATIONS

Image

SKU

Price

**Stock
Status** **Stock
Quantity**

Description

Size

FDA Class II Medical Device



PR515-MLA \$3,439.00

Experience true weightlessness. The Cloud™ is a luxury power lift recliner offering our patented Twilight Tilt Technology & MaxiComfort® Positioning for long-term comfort and support.

Med/Large

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
FDA Class II Medical Device 	PR515	\$3,439.00			Experience true weightlessness. The Cloud™ is a luxury power lift recliner offering our patented Twilight Tilt Technology & MaxiComfort® Positioning for long-term comfort and support.	Small/Med

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Med/Large, Small/Med
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	Twilight Tilt
Comfort Zone(s)	5
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Quick Release
Head Rest Cover	Standard
Arm Rest Covers	No
Overall Width	35.5"



Overall Height	44"
Seat Width	23.5"
Floor to Top of Seat	20"
Seat Depth	18.5"
Seat to Top of Back	28"
Seat to Top of Arm	8"
Distance Required from Wall Reclined	23"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	Yes
Nirvana Air Massage & Heat Available	No
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	Yes

Straight Lift Option	Standard
USB Charging Port	Control Box
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA