

FDA Class II Medical Device



CLOUD CLASSIC

\$2,479.00 - \$2,839.00 Price range:
\$2,479.00 through \$2,839.00

The Cloud™ Classic is an *extremely* comfortable and popular power lift recliner offering our MaxiComfort® Positioning & custom bucket seat design for ultimate comfort and support.

SKU: PR510

Categories: [2 Zone Chairs](#), [Heavy Duty](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

VARIATIONS

FDA Class II Medical Device



PR510-MLA \$2,479.00


Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
	PR510-MLA	\$2,479.00			Relax and recharge in comfort and style with the PR517 Cloud Power Lift Recliner - 4-Zone Comfort Zone positioning system!	Med/Large

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
<small>FDA Class II Medical Device</small> 	PR510-MXW	\$2,839.00			<p>The Cloud™ is a popular Golden power lift recliner offering our MaxiComfort® Positioning & 600 lbs. weight capacity.</p> <p>The Cloud™ is an extremely comfortable and popular power lift recliner offering our MaxiComfort® Positioning & custom bucket seat design for ultimate comfort and support.</p>	Medium/Wide
<small>FDA Class II Medical Device</small> 	PR510-SME	\$2,479.00			<p>The Cloud™ is an extremely comfortable and popular power lift recliner offering our MaxiComfort® Positioning & custom bucket seat design for ultimate comfort and support.</p>	Small/Med

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Med/Large, Medium/Wide, Small/Med
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No

Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Biscuit
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	No
Overall Width	33.5"
Overall Height	41"
Seat Width	18.5"
Floor to Top of Seat	19"
Seat Depth	19.5"
Seat to Top of Back	27"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	23"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes

HeatWave Technology Option	Yes
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	No
Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA