

FDA Class II Medical Device



COMFORTER LARGE WITH ZG+

\$2,959.00

Available in three sizes, this Comforter power lift recliner offers our patented Twilight Tilt Technology, MaxiComfort® and Zero Gravity+ Positioning!

SKU: PR545-LAR

Categories: [3 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Small to Tall](#), [Twilight](#)

ADDITIONAL INFORMATION

Size	Large
Class II Medical Device	YES
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	3
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	Standard



Arm Rest Covers	Standard
Overall Width	34"
Overall Height	43
Seat Width	23"
Floor to Top of Seat	21"
Seat Depth	21"
Seat to Top of Back	30"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	24"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	NO
Foot Rest Extension Option	Yes



Straight Lift Option	Standard
USB Charging Port	Control Box
Battery Backup	Lithium-Ion
SmartTek® Diagnostic System	Standard
Place of Origin	USA