



COMFORTER LARGE WITH ZG+

\$2,839.00

Available in three sizes, this Comforter power lift recliner offers our patented Twilight Tilt Technology, MaxiComfort® and Zero Gravity+ Positioning!

SKU: PR545-LAR

Categories: [3 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Small to Tall](#), [Twilight](#)

Tags: [New](#)

ADDITIONAL INFORMATION

| | |
|---------------------------------|-------------------------|
| Size | Large |
| Class II Medical Device | YES |
| Weight Capacity | 375 lbs |
| Positioning | MaxiComfort® |
| Comfort Zone(s) | 3 |
| Twilight Tilt Technology | Yes |
| ZG+ | Yes |
| Chaise Seat | Yes |
| Pocketed Coil Seat | Yes |
| Back Type | Seam |
| Removable Back Type | Quick Release |
| Head Rest Cover | Standard |

| | |
|---|---------------|
| Arm Rest Covers | Standard |
| Overall Width | 34" |
| Overall Height | 43 |
| Seat Width | 23" |
| Floor to Top of Seat | 21" |
| Seat Depth | 21" |
| Seat to Top of Back | 30" |
| Seat to Top of Arm | 6" |
| Distance Required from Wall Reclined | 24" |
| Extension Reclined from Back of Seat to End of Foot Rest | 43" |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | Yes |
| Warranty on Electrical Parts | 3 years |
| Pro Rated Years 4 Thru 7 | Yes |
| Warranty on Mechanical Labor | 3 years |
| Warranty on Electrical Labor | 1 year |
| Hand Control Routing | Inside Pocket |
| Extra Pocket Option | Standard |
| Heat & Vibrating Massage Option | Yes |
| HeatWave Technology Option | No |
| Articulating Headrest Option | NO |
| Power Adjustable Lumbar Option | NO |
| Foot Rest Extension Option | Yes |

| | |
|------------------------------------|-------------|
| Straight Lift Option | Standard |
| USB Charging Port | Control Box |
| Battery Backup | Lithium-Ion |
| SmartTek® Diagnostic System | Standard |
| Place of Origin | USA |