



COMFORTER LARGE

\$1,599.00

Whether you're 5'0" to 6'8" and 300 to 375 lbs., the Comforter Series from Golden Technologies has a power lift chair in just your size—actually, five sizes. Featuring a plush seamed backrest with lumbar support, a soft, comfortable seat and our "grip rite" arm design, Comforter lift recliners provide the perfect fit.

SKU: PR531-LAR

Categories: [1 Zone Chairs](#), [Comforter](#), [Heavy Duty](#), [Lift Recliners](#), [Small to Tall](#)

ADDITIONAL INFORMATION

Size	Large
Class II Medial Device	Yes
Weight Capacity	375 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Bolt

Head Rest Cover	Standard
Arm Rest Covers	Standard
Overall Width	34"
Overall Height	44"
Seat Width	23"
Floor to Top of Seat	21"
Seat Depth	20.5"
Seat to Top of Back	28.5"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	19"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No

Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	Hand Control
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA