



COMFORTER MEDIUM WIDE

\$2,049.00

Golden Technologies recognizes that “One size does not fit all” when it comes to power lift chairs. We’ve expanded our classic Comforter Series with wider seating areas and greater lifting capacities. This medium-wide power lift recliner has a weight capacity of 500 pounds. Most of the five lift recliners in this series use two heavy-duty motors for strength and support, with a three-position style recline mechanism.

SKU: PR531-MXW

Categories: [1 Zone Chairs](#), [Comfort Zones](#), [Comforter Wides](#), [Heavy Duty](#), [Lift Recliners](#), [Small to Tall](#)

ADDITIONAL INFORMATION

Size	Medium/Wide
Class II Medial Device	Yes
Weight Capacity	500 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes

Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	Standard
Overall Width	38"
Overall Height	43.5"
Seat Width	27"
Floor to Top of Seat	20.5"
Seat Depth	20.5"
Seat to Top of Back	29"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	18"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No

Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	Hand Control
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA