



**COMFORTER MEDIUM WIDE**

**\$2,049.00**

Golden Technologies recognizes that “One size does not fit all” when it comes to power lift chairs. We’ve expanded our classic Comforter Series with wider seating areas and greater lifting capacities. This medium-wide power lift recliner has a weight capacity of 500 pounds. Most of the five lift recliners in this series use two heavy-duty motors for strength and support, with a three-position style recline mechanism.

---

**SKU:** PR531-MXW

**Categories:** [1 Zone Chairs](#), [3-Position](#), [Comfort Zones](#), [Comforter Wides](#), [Lift Recliners](#), [Small to Tall](#)

**ADDITIONAL INFORMATION**

<b>Size</b>	Medium/Wide
<b>Class II Medial Device</b>	Yes
<b>Weight Capacity</b>	<a href="#">500 lbs</a>
<b>Positioning</b>	3 Position
<b>Comfort Zone(s)</b>	1
<b>Twilight Tilt Technology</b>	No
<b>ZG+</b>	No
<b>Chaise Seat</b>	Yes

<b>Pocketed Coil Seat</b>	No
<b>Back Type</b>	Seam
<b>Removable Back Type</b>	Bolt
<b>Head Rest Cover</b>	Standard
<b>Arm Rest Covers</b>	Standard
<b>Overall Width</b>	38"
<b>Overall Height</b>	43.5"
<b>Seat Width</b>	27"
<b>Floor to Top of Seat</b>	20.5"
<b>Seat Depth</b>	20.5"
<b>Seat to Top of Back</b>	29"
<b>Seat to Top of Arm</b>	7.5"
<b>Distance Required from Wall Reclined</b>	18"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	42"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty on Electrical Parts</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	No

<b>Articulating Headrest Option</b>	NO
<b>Power Adjustable Lumbar Option</b>	No
<b>Foot Rest Extension Option</b>	Yes
<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	Hand Control
<b>Battery Backup</b>	9 volt
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA