



COMFORTER MEDIUM WIDE

\$2,049.00

Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.

SKU: PR531-MXW

Categories: [1 Zone Chairs](#), [Comfort Zones](#), [Comforter Wides](#), [Heavy Duty](#), [Lift Recliners](#), [Small to Tall](#)

ADDITIONAL INFORMATION

| | |
|---------------------------------|-------------------------|
| Size | Medium/Wide |
| Class II Medical Device | Yes |
| User Weight Capacity | 500 lbs |
| Positioning | 3 Position |
| Comfort Zone(s) | 1 |
| Twilight Tilt Technology | No |
| ZG+ | No |
| Chaise Seat | Yes |
| Pocketed Coil Seat | No |
| Back Type | Seam |
| Removable Back Type | Bolt |
| Head Rest Cover | Standard |
| Arm Rest Covers | Standard |



| | |
|---|---------------|
| Overall Width | 38" |
| Overall Height | 43.5" |
| Seat Width | 27" |
| Floor to Top of Seat | 20.5" |
| Seat Depth | 20.5" |
| Seat to Top of Back | 29" |
| Seat to Top of Arm | 7.5" |
| Distance Required from Wall Reclined | 18" |
| Extension Reclined from Back of Seat to End of Foot Rest | 42" |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | Yes |
| Warranty Electronics | 3 years |
| Pro Rated Years 4 Thru 7 | Yes |
| Warranty on Mechanical Labor | 3 years |
| Warranty on Electrical Labor | 1 year |
| Hand Control Routing | Inside Pocket |
| Extra Pocket Option | Standard |
| Heat & Vibrating Massage Option | Yes |
| HeatWave Technology Option | No |
| Articulating Headrest Option | NO |
| Power Adjustable Lumbar Option | No |
| Foot Rest Extension Option | Yes |
| Straight Lift Option | Standard |



| | |
|-----------------------------------|--------------|
| USB Charging Port | Hand Control |
| Battery Backup | 9 volt |
| SmartTek Diagnostic System | Standard |
| Place of Origin | USA |