



COMFORTER MEDIUM

\$1,449.00

Golden Technologies recognizes that “One size does not fit all” when it comes to power lift recliners. Featuring a plush seamed backrest with lumbar support, a soft, comfortable seat and our “grip rite” arm design, the Comforter power lift recliner gives you the perfect fit. This stylish medium power lift recliner has a weight capacity of 375 pounds and is available in numerous other sizes.

SKU: PR531-MED

Categories: [1 Zone Chairs](#), [Comforter](#), [Heavy Duty](#), [Lift Recliners](#), [Small to Tall](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medial Device	Yes
Weight Capacity	375 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam

Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	Standard
Overall Width	32"
Overall Height	41"
Seat Width	21"
Floor to Top of Seat	20"
Seat Depth	19"
Seat to Top of Back	27"
Seat to Top of Arm	5"
Distance Required from Wall Reclined	17"
Extension Reclined from Back of Seat to End of Foot Rest	41"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO

Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	Hand Control
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA