

FDA Class II Medical Device



COMFORTER PETITE SMALL

\$1,499.00

Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.

SKU: PR531-PSA

Categories: [1 Zone Chairs](#), [Comforter](#), [Heavy Duty](#), [Lift Recliners](#), [Small to Tall](#)

ADDITIONAL INFORMATION

Size	Petite/Small
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	Standard



Overall Width	31.5"
Overall Height	40"
Seat Width	21"
Floor to Top of Seat	19"
Seat Depth	19"
Seat to Top of Back	27"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	16"
Extension Reclined from Back of Seat to End of Foot Rest	41"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard

USB Charging Port	Hand Control
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA