



**COMFORTER PETITE SMALL**

**\$1,449.00**

Whether you're 5'0" to 6'8" and 300 to 375 lbs., the Comforter Series lift recliners from Golden Technologies has just your size—actually, five sizes. Featuring a plush seamed backrest with lumbar support, a soft, comfortable seat and our "grip rite" arm design, the Comforter power lift recliner gives you the perfect fit.

---

**SKU:** PR531-PSA

**Categories:** [1 Zone Chairs](#), [Comforter](#), [Heavy Duty](#), [Lift Recliners](#), [Small to Tall](#)

**ADDITIONAL INFORMATION**

<b>Size</b>	Petite/Small
<b>Class II Medial Device</b>	Yes
<b>Weight Capacity</b>	<a href="#">375 lbs</a>
<b>Positioning</b>	3 Position
<b>Comfort Zone(s)</b>	1
<b>Twilight Tilt Technology</b>	No
<b>ZG+</b>	No
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	Yes
<b>Back Type</b>	Seam

<b>Removable Back Type</b>	Bolt
<b>Head Rest Cover</b>	Standard
<b>Arm Rest Covers</b>	Standard
<b>Overall Width</b>	31.5"
<b>Overall Height</b>	40"
<b>Seat Width</b>	21"
<b>Floor to Top of Seat</b>	19"
<b>Seat Depth</b>	19"
<b>Seat to Top of Back</b>	27"
<b>Seat to Top of Arm</b>	6"
<b>Distance Required from Wall Reclined</b>	16"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	41"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty on Electrical Parts</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	No
<b>Articulating Headrest Option</b>	NO

<b>Power Adjustable Lumbar Option</b>	No
<b>Foot Rest Extension Option</b>	Yes
<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	Hand Control
<b>Battery Backup</b>	9 volt
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA