



COMFORTER PETITE SMALL

\$1,498.00

Whether you're 5'0" to 6'8" and 300 to 375 lbs., the Comforter Series lift recliners from Golden Technologies has just your size—actually, five sizes. Featuring a plush seamed backrest with lumbar support, a soft, comfortable seat and our "grip rite" arm design, the Comforter power lift recliner gives you the perfect fit.

SKU: PR531-PSA

Categories: [1 Zone Chairs](#), [3-Position, Comforter](#), [Lift Recliners](#), [Small to Tall](#)

ADDITIONAL INFORMATION

Size	Petite/Small
Weight capacity	375
Positioning	3 Position
Twilight Tilt	No
Chaise seat	Yes
Pocketed coil seat	Yes
Back type	Seam
Removable Back Type	Bolt
Head rest cover	YES
Arm rest covers	Yes

Overall width	32"
Overall height	40"
Width between arms	21"
Floor to top of seat	20"
Seat depth	19"
Seat to top of back	27"
Seat to top of arm	6"
Dis req from wall reclined	16"
Ext reclined from back	41"
Life time warranty frame	Yes
Warranty electrical parts	3 Yrs.
Pro rated 4 thru 7	Yes
Warranty mechanical labor	3 Yrs.
Warranty electrical labor	1 Yr.
Hand control routing	Inside Arm
Extra pocket option	Yes
Heat massage option	Yes
Adjustable headrest option	No
Power adjustable lumbar	No
Foot rest ext option	Yes
Straight lift option	Yes
Heatwave technology option	No
Usb charging port	YES
Place of origin	USA