



## COMFORTER TALL WITH ZG+

**\$3,089.00**

Available in three sizes, this Comforter power lift recliner offers our patented Twilight Tilt Technology, MaxiComfort® and Zero Gravity+ Positioning!

---

**SKU:** PR545-TAL

**Categories:** [3 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Small to Tall](#), [Twilight](#)

**Tags:** [New](#)

### ADDITIONAL INFORMATION

<b>Size</b>	Tall
<b>Class II Medical Device</b>	YES
<b>Weight Capacity</b>	<a href="#">375 lbs</a>
<b>Positioning</b>	MaxiComfort®
<b>Comfort Zone(s)</b>	3
<b>Twilight Tilt Technology</b>	Yes
<b>ZG+</b>	Yes
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	Yes
<b>Back Type</b>	Seam
<b>Removable Back Type</b>	Quick Release
<b>Head Rest Cover</b>	Standard

<b>Arm Rest Covers</b>	Standard
<b>Overall Width</b>	34"
<b>Overall Height</b>	48"
<b>Seat Width</b>	23"
<b>Floor to Top of Seat</b>	23"
<b>Seat Depth</b>	22"
<b>Seat to Top of Back</b>	32"
<b>Seat to Top of Arm</b>	6"
<b>Distance Required from Wall Reclined</b>	28"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	45"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty on Electrical Parts</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	No
<b>Articulating Headrest Option</b>	NO
<b>Power Adjustable Lumbar Option</b>	NO
<b>Foot Rest Extension Option</b>	Yes

<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	Control Box
<b>Battery Backup</b>	Lithium-Ion
<b>SmartTek® Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA