

FDA Class II Medical Device



COMFORTER

\$1,499.00 - \$2,639.00 Price range:
\$1,499.00 through \$2,639.00

Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.

SKU: PR531

Categories: [1 Zone Chairs](#), [Comfort Zones](#), [Comforter Wides](#), [Heavy Duty](#), [Lift Recliners](#), [Small to Tall](#)

VARIATIONS

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
FDA Class II Medical Device 	PR531-LAR	\$1,649.00			Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.	Large
FDA Class II Medical Device 	PR531-MED	\$1,499.00			Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.	Medium

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
FDA Class II Medical Device 	PR531-MXW	\$2,119.00			Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.	Medium/Wide
FDA Class II Medical Device 	PR531-PSA	\$1,499.00			Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.	Petite/Small
FDA Class II Medical Device 	PR531-TAL	\$1,759.00			Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.	Tall
FDA Class II Medical Device 	PR531-TXW	\$2,639.00			Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.	Tall/Wide
FDA Class II Medical Device 	PR531-SXW	\$1,819.00			Small to Tall, We Build Them All! Golden's Comforter Series of lift recliners has just your size.	Small/Wide

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Large, Medium, Medium/Wide, Petite/Small, Small/Wide, Tall, Tall/Wide
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	Standard
Overall Width	34"
Overall Height	40.5"
Seat Width	23"
Floor to Top of Seat	19.5"
Seat Depth	20"
Seat to Top of Back	27"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	16"
Extension Reclined from Back of Seat to End of Foot Rest	40"

Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	Hand Control
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA