



DAY DREAMER MEDIUM

\$2,899.00

Dream away in the DayDreamer™ by Golden. This chair offers our MaxiComfort® Positioning and an exclusive PowerPillow™ articulating headrest and power lumbar.

SKU: PR632-MED

Categories: [4 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medial Device	Yes
Weight Capacity	375 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	4
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Bolt
Head Rest Cover	No
Arm Rest Covers	No

Overall Width	38"
Overall Height	45"
Seat Width	22.5"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	31"
Seat to Top of Arm	6.5"
Distance Required from Wall Reclined	22"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	No
HeatWave Technology Option	No
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	Yes
Straight Lift Option	Standard

USB Charging Port

Control Box

Battery Backup

Lithium-ion

SmartTek Diagnostic System

Standard

Place of Origin

USA