



ELARA LARGE

\$1,239.00

Inspired by the moon, the DeLuna™ Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep.

SKU: PR118-LAR

Categories: [1 Zone Chairs](#), [DeLuna](#), [Heavy Duty](#), [Lift Recliners](#)

ADDITIONAL INFORMATION

Size	Large
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	3 Position
Comfort Zone(s)	1
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	Standard
Arm Rest Covers	No



Overall Width	34"
Overall Height	46"
Seat Width	23"
Floor to Top of Seat	19.5"
Seat Depth	21"
Seat to Top of Back	31"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	20"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Arm
Extra Pocket Option	No
Heat & Vibrating Massage Option	No
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	No
Straight Lift Option	Standard

USB Charging Port	Hand Control
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	Imported