

FDA Class II Medical Device



EZ SLEEPER SLIM WITH TWILIGHT

\$3,259.00

The 5-Zone EZ Sleeper Slim features Golden's patented Twilight Tilt Technology and MaxiComfort Positioning. Upgrade with Nirvana Swedish Air Massage & FIR Heat!

SKU: PR764-MED

Categories: [5 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No



Arm Rest Covers	No
Overall Width	33.5"
Overall Height	45.5"
Seat Width	22.5"
Floor to Top of Seat	20"
Seat Depth	22"
Seat to Top of Back	30"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	21.5"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
Nirvana Air Massage & Heat Available	Yes
HeatWave Technology Option	Yes
Articulating Headrest Option	Standard



Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	No
Straight Lift Option	Standard
USB Charging Port	YES
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA