

FDA Class II Medical Device



## EZ SLEEPER SLIM

The 4-Zone EZ Sleeper Slim has a modern look with a tailored back and features Golden's patented MaxiComfort Positioning.

---

**SKU:** PR750-MED

**Categories:** [4 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

**Tags:** [New](#)

## ADDITIONAL INFORMATION

<b>Size</b>	Medium
<b>Class II Medical Device</b>	Yes
<b>User Weight Capacity</b>	<a href="#">400 lbs</a>
<b>Positioning</b>	MaxiComfort®
<b>Twilight Tilt Technology</b>	No
<b>ZG+</b>	Yes
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	Yes
<b>Back Type</b>	Seam
<b>Removable Back Type</b>	Quick Release
<b>Head Rest Cover</b>	No
<b>Arm Rest Covers</b>	No
<b>Overall Width</b>	33.5"



<b>Overall Height</b>	45.5"
<b>Seat Width</b>	22.5"
<b>Floor to Top of Seat</b>	20"
<b>Seat Depth</b>	22"
<b>Seat to Top of Back</b>	30"
<b>Seat to Top of Arm</b>	7"
<b>Distance Required from Wall Reclined</b>	21.5"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	42"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty Electronics</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>Nirvana Air Massage &amp; Heat Available</b>	Yes
<b>HeatWave Technology Option</b>	Yes
<b>Articulating Headrest Option</b>	Standard
<b>Power Adjustable Lumbar Option</b>	Standard
<b>Foot Rest Extension Option</b>	No



<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	YES
<b>Battery Backup</b>	Lithium-ion
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA