

FDA Class II Medical Device



EZ SLEEPER WITH TWILIGHT MEDIUM/LARGE

\$3,299.00

Get the rest and comfort you need with The EZ Sleeper! This modern power lift recliner offers our patented Twilight Tilt Technology & MaxiComfort® Positioning.

SKU: PR761-MLA

Categories: [5 Zone Chairs](#), [Lift Recliners](#), [Patented MaxiComfort Positioning](#), [Patented Twilight Tilt](#), [Twilight](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Size	Med/Large
Class II Medical Device	Yes
Weight Capacity	400 lbs
Positioning	Twilight Tilt
Comfort Zone(s)	5
Twilight Tilt Technology	Yes
ZG+	Yes
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	No



Arm Rest Covers	No
Overall Width	34"
Overall Height	46"
Seat Width	22"
Floor to Top of Seat	21"
Seat Depth	21"
Seat to Top of Back	30"
Seat to Top of Arm	8"
Distance Required from Wall Reclined	24"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	Yes
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	Yes



Straight Lift Option	Standard
USB Charging Port	Control Box
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA