



EZ SLEEPER MEDIUM LARGE

\$2,499.00

The 2-Zone EZ Sleeper Lift Chair is just for you! This luxury power lift recliner combines a contemporary design with our patented MaxiComfort positioning for unmatched style and comfort. Its ultramodern, clean lines pair perfectly with contemporary décor, especially in smaller spaces.

SKU: PR735-MLA

Categories: [2 Zone Chairs](#), [Comfort Zones](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

ADDITIONAL INFORMATION

Size	Med/Large
Class II Medial Device	Yes
Weight Capacity	375 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam

Removable Back Type	Bolt
Head Rest Cover	No
Arm Rest Covers	No
Overall Width	34.5"
Overall Height	45.5"
Seat Width	21.5"
Floor to Top of Seat	20"
Seat Depth	22.5"
Seat to Top of Back	30"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	20"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO

Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA