



**COMFORTER WITH MAXICOMFORT
MEDIUM/WIDE**

\$2,749.00

Small to Tall, We Build Them All! Find your perfect fit with the Medium Wide Comforter featuring Golden's patented and popular MaxiComfort® Positioning.

SKU: PR535-MXW

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Small to Tall](#)

ADDITIONAL INFORMATION

Size	Medium/Wide
Class II Medial Device	Yes
Weight Capacity	500 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	Standard

Arm Rest Covers	Standard
Overall Width	38"
Overall Height	45"
Seat Width	27"
Floor to Top of Seat	20"
Seat Depth	20"
Seat to Top of Back	30"
Seat to Top of Arm	5"
Distance Required from Wall Reclined	25"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes

Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA