

FDA Class II Medical Device



COMFORTER WITH MAXICOMFORT PETITE SMALL

\$2,199.00

Small to Tall, We Build Them All! Find your perfect fit with the Petite Small Comforter featuring Golden's popular MaxiComfort® Positioning.

SKU: PR535-PSA

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Small to Tall](#)

Tags: [Customer Favorite](#)

ADDITIONAL INFORMATION

Size	Petite/Small
Class II Medical Device	Yes
Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	Standard

Arm Rest Covers	Standard
Overall Width	32"
Overall Height	40"
Seat Width	21"
Floor to Top of Seat	20"
Seat Depth	19"
Seat to Top of Back	26"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	22"
Extension Reclined from Back of Seat to End of Foot Rest	40"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Foot Rest Extension Option	Yes
Straight Lift Option	Standard

USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA