



COMFORTER WITH MAXICOMFORT TALL

\$2,579.00

Small to Tall, We Build Them All! Find your perfect fit with the Tall Comforter featuring Golden's patented and popular MaxiComfort® Positioning.

SKU: PR535-TAL

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Small to Tall](#)

Tags: [Customer Favorite](#)

ADDITIONAL INFORMATION

Size	Tall
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	Standard



Arm Rest Covers	Standard
Overall Width	34"
Overall Height	48"
Seat Width	23"
Floor to Top of Seat	23"
Seat Depth	21"
Seat to Top of Back	31"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	26"
Extension Reclined from Back of Seat to End of Foot Rest	45"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes



Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA