

FDA Class II Medical Device



COMFORTER WITH MAXICOMFORT

\$2,269.00 - \$2,989.00 Price range:
\$2,269.00 through \$2,989.00

Small to Tall, We Build Them All! Find your perfect fit with the Petite Small Comforter featuring Golden's popular MaxiComfort® Positioning.

SKU: PR535

Categories: [2 Zone Chairs](#), [Heavy Duty](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#), [Small to Tall](#)

Tags: [Customer Favorite](#)

VARIATIONS

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
	PR535-LAR	\$2,319.00			Small to Tall, We Build Them All! Find your perfect fit with the Large Comforter featuring Golden's patented and popular MaxiComfort® Positioning.	Large

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
<p><small>FDA Class II Medical Device</small></p> 	PR535-LXW	\$2,989.00			<p>Small to Tall, We Build Them All! Now with a 600 lb. weight capacity. Find your perfect fit with the Large Wide Comforter featuring Golden's patented and popular MaxiComfort® Positioning.</p>	Large/Wide
	PR535-MED	\$2,269.00			<p>Small to Tall, We Build Them All! Find your perfect fit with the Medium Comforter featuring Golden's popular MaxiComfort® Positioning.</p>	Medium
	PR535-MXW	\$2,839.00			<p>Small to Tall, We Build Them All! Find your perfect fit with the Medium Wide Comforter featuring Golden's patented and popular MaxiComfort® Positioning.</p>	Medium/Wide

Image	SKU	Price	Stock Status	Stock Quantity	Description	Size
	PR535-TAL	\$2,579.00			<p>Small to Tall, We Build Them All! Find your perfect fit with the Tall Comforter featuring Golden's patented and popular MaxiComfort® Positioning.</p>	Tall
<small>FDA Class II Medical Device</small> 	PR535-PSA	\$2,269.00			<p>Small to Tall, We Build Them All! Find your perfect fit with the Petite Small Comforter featuring Golden's popular MaxiComfort® Positioning.</p>	Petite/Small

ADDITIONAL INFORMATION

Dimensions	N/A
Size	Large, Large/Wide, Medium, Medium/Wide, Petite/Small, Tall
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No

ZG+	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	Standard
Arm Rest Covers	Standard
Overall Width	32"
Overall Height	40"
Seat Width	21"
Floor to Top of Seat	20"
Seat Depth	19"
Seat to Top of Back	26"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	22"
Extension Reclined from Back of Seat to End of Foot Rest	40"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes

HeatWave Technology Option	No
Articulating Headrest Option	NO
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA