



**REGAL - REPLACED BY CLOUD+**

**\$2,899.00**

Explore the NEW Cloud+ power lift recliner - Available in [Small](#), [Medium](#) & [Large](#).

---

**SKU:** PR504-MLA-ST4

**Categories:** [4 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

**ADDITIONAL INFORMATION**

<b>Size</b>	Med/Large
<b>Class II Medial Device</b>	Yes
<b>Weight Capacity</b>	<a href="#">400 lbs</a>
<b>Positioning</b>	MaxiComfort®
<b>Comfort Zone(s)</b>	4
<b>Twilight Tilt Technology</b>	No
<b>ZG+</b>	No
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	No
<b>Back Type</b>	Seam
<b>Removable Back Type</b>	Quick Release
<b>Head Rest Cover</b>	No
<b>Arm Rest Covers</b>	No

<b>Overall Width</b>	40"
<b>Overall Height</b>	48"
<b>Seat Width</b>	23
<b>Floor to Top of Seat</b>	20.5"
<b>Seat Depth</b>	21.5"
<b>Seat to Top of Back</b>	33"
<b>Seat to Top of Arm</b>	7.5"
<b>Distance Required from Wall Reclined</b>	22"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	42.5"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty on Electrical Parts</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Right Arm Compartment
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	Standard
<b>Articulating Headrest Option</b>	Standard
<b>Power Adjustable Lumbar Option</b>	Standard
<b>Foot Rest Extension Option</b>	Yes
<b>Straight Lift Option</b>	Standard

<b>USB Charging Port</b>	Left Arm Compartment
<b>Battery Backup</b>	Lithium-ion
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA