



REGAL MEDIUM LARGE

\$2,899.00

Fit for the King or Queen of the House! The Regal has MaxiComfort® Positioning, a higher back, numerous accessories & comes standard with HeatWave™ Technology.

SKU: PR504-MLA-ST4

Categories: [4 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Zero Gravity \(MaxiComfort\)](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Size	Med/Large
Class II Medial Device	YES
Weight Capacity	400
Positioning	MaxiComfort
Comfort Zone(s)	4
Twilight Tilt Technology	NO
ZG+	NO
Chaise Seat	YES
Pocketed Coil Seat	NO
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	NO
Arm Rest Covers	No



Overall Width	40"
Overall Height	48"
Seat Width	23
Floor to Top of Seat	20.5"
Seat Depth	21.5"
Seat to Top of Back	33"
Seat to Top of Arm	7.5"
Distance Required from Wall Reclined	22"
Extension Reclined from Back of Seat to End of Foot Rest	42.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Right Arm Compartment
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	YES
HeatWave Technology Option	Standard
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	YES
Straight Lift Option	Standard



USB Charging Port	Left Arm Compartment
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	USA