

FDA Class II Medical Device



RELAXER LARGE

\$2,499.00

The Relaxer features an adjustable three-pillow waterfall back & MaxiComfort Positioning. Available in three sizes.

SKU: PR766-LAR

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

ADDITIONAL INFORMATION

Size	Large
Class II Medical Device	Yes
User Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	3 Pillow
Removable Back Type	Bolt
Head Rest Cover	Standard
Arm Rest Covers	Standard
Overall Width	38.5"

Overall Height	48"
Seat Width	22.5"
Floor to Top of Seat	22"
Seat Depth	22.5"
Seat to Top of Back	29.5"
Seat to Top of Arm	5.5"
Distance Required from Wall Reclined	27"
Extension Reclined from Back of Seat to End of Foot Rest	47"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty Electronics	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	NO

Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA