



RELAXER LARGE

\$2,649.00

The “Relaxer” is Golden’s original MaxiComfort power lift chair recliner offering a three pillow waterfall back and a deep seat for lateral support. This comfy recliner is available in two sizes and each of the three back pillows can be customized by the user to find that perfect fit.

SKU: PR766-LAR

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Zero Gravity \(MaxiComfort\)](#)

ADDITIONAL INFORMATION

Size	Large
Class II Medial Device	YES
Weight Capacity	375
Positioning	MaxiComfort
Comfort Zone(s)	2
Twilight Tilt Technology	NO
Chaise Seat	YES
Pocketed Coil Seat	YES
Back Type	3 Pillow
Removable Back Type	Bolt
Head Rest Cover	Standard



Arm Rest Covers	Standard
Overall Width	38.5"
Overall Height	48"
Seat Width	22.5"
Floor to Top of Seat	22"
Seat Depth	22.5"
Seat to Top of Back	29.5"
Seat to Top of Arm	5.5"
Distance Required from Wall Reclined	27"
Extension Reclined from Back of Seat to End of Foot Rest	47"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	YES
HeatWave Technology Option	NO
Articulating Headrest Option	NO
Power Adjustable Lumbar Option	NO
Foot Rest Extension Option	YES



Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA