



## RELAXER LARGE

**\$2,499.00**

The “Relaxer” is Golden’s original MaxiComfort power lift chair recliner offering a three pillow waterfall back and a deep seat for lateral support. This comfy recliner is available in two sizes and each of the three back pillows can be customized by the user to find that perfect fit.

---

**SKU:** PR766-LAR

**Categories:** [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

### ADDITIONAL INFORMATION

<b>Size</b>	Large
<b>Class II Medial Device</b>	Yes
<b>Weight Capacity</b>	<a href="#">375 lbs</a>
<b>Positioning</b>	MaxiComfort®
<b>Comfort Zone(s)</b>	2
<b>Twilight Tilt Technology</b>	No
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	Yes
<b>Back Type</b>	3 Pillow
<b>Removable Back Type</b>	Bolt

<b>Head Rest Cover</b>	Standard
<b>Arm Rest Covers</b>	Standard
<b>Overall Width</b>	38.5"
<b>Overall Height</b>	48"
<b>Seat Width</b>	22.5"
<b>Floor to Top of Seat</b>	22"
<b>Seat Depth</b>	22.5"
<b>Seat to Top of Back</b>	29.5"
<b>Seat to Top of Arm</b>	5.5"
<b>Distance Required from Wall Reclined</b>	27"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	47"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty on Electrical Parts</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	No
<b>Articulating Headrest Option</b>	NO
<b>Power Adjustable Lumbar Option</b>	No

<b>Foot Rest Extension Option</b>	Yes
<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	NO
<b>Battery Backup</b>	9 volt
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA