



## RELAXER MEDIUM

**\$2,399.00**

The Relaxer is the original power lift chair recliner featuring our patented MaxiComfort® Ultimate Recline Technology standard! The Relaxer lift chair is a true Golden customer favorite, now with a design refresh with modern, well-appointed features like super plush rounded armrests, an adjustable three-pillow waterfall back, an updated pocket style, and an extra-wide chaise pad for better lower leg support.

---

**SKU:** PR766-MED

**Categories:** [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

## ADDITIONAL INFORMATION

<b>Size</b>	Medium
<b>Class II Medial Device</b>	Yes
<b>Weight Capacity</b>	<a href="#">375 lbs</a>
<b>Positioning</b>	MaxiComfort®
<b>Comfort Zone(s)</b>	2
<b>Twilight Tilt Technology</b>	No
<b>Chaise Seat</b>	Yes
<b>Pocketed Coil Seat</b>	Yes

<b>Back Type</b>	3 Pillow
<b>Removable Back Type</b>	Bolt
<b>Head Rest Cover</b>	Yes
<b>Arm Rest Covers</b>	Standard
<b>Overall Width</b>	36"
<b>Overall Height</b>	44.5"
<b>Seat Width</b>	22.5"
<b>Floor to Top of Seat</b>	20"
<b>Seat Depth</b>	20.5"
<b>Seat to Top of Back</b>	28.5"
<b>Seat to Top of Arm</b>	5"
<b>Distance Required from Wall Reclined</b>	26"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	44"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	Yes
<b>Warranty on Electrical Parts</b>	3 years
<b>Pro Rated Years 4 Thru 7</b>	Yes
<b>Warranty on Mechanical Labor</b>	3 years
<b>Warranty on Electrical Labor</b>	1 year
<b>Hand Control Routing</b>	Inside Pocket
<b>Extra Pocket Option</b>	Standard
<b>Heat &amp; Vibrating Massage Option</b>	Yes
<b>HeatWave Technology Option</b>	No

<b>Articulating Headrest Option</b>	NO
<b>Power Adjustable Lumbar Option</b>	No
<b>Foot Rest Extension Option</b>	Yes
<b>Straight Lift Option</b>	Standard
<b>USB Charging Port</b>	NO
<b>Battery Backup</b>	9 volt
<b>SmartTek Diagnostic System</b>	Standard
<b>Place of Origin</b>	USA