

FDA Class II Medical Device



RELAXER MEDIUM

\$2,399.00

The Relaxer is the original power lift chair recliner featuring our patented MaxiComfort® Ultimate Recline Technology standard! The Relaxer lift chair is a true Golden customer favorite, now with a design refresh with modern, well-appointed features like super plush rounded armrests, an adjustable three-pillow waterfall back, an updated pocket style, and an extra-wide chaise pad for better lower leg support.

SKU: PR766-MED

Categories: [2 Zone Chairs](#), [Lift Recliners](#), [MaxiComfort](#), [Patented MaxiComfort Positioning](#)

ADDITIONAL INFORMATION

Size	Medium
Class II Medical Device	Yes
Weight Capacity	400 lbs
Positioning	MaxiComfort®
Comfort Zone(s)	2
Twilight Tilt Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	Yes

Back Type	3 Pillow
Removable Back Type	Bolt
Head Rest Cover	Yes
Arm Rest Covers	Standard
Overall Width	36"
Overall Height	44.5"
Seat Width	22.5"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	28.5"
Seat to Top of Arm	5"
Distance Required from Wall Reclined	26"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 years
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 years
Warranty on Electrical Labor	1 year
Hand Control Routing	Inside Pocket
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	Yes
HeatWave Technology Option	No

Articulating Headrest Option	NO
Power Adjustable Lumbar Option	No
Foot Rest Extension Option	Yes
Straight Lift Option	Standard
USB Charging Port	NO
Battery Backup	9 volt
SmartTek Diagnostic System	Standard
Place of Origin	USA